<div style="display: flex; justify-content: space-between;"> 5293140602 Adhesive Resin Cement </div> <div style="display: flex; justify-content: space-between;"> Komposit-Befestigungsmaterial für adhäsive Verankerungen Résine de collage </div> <div style="display: flex; justify-content: space-between;"> Cemento resina adhesivo Cemento Adhesivo para Materiales de Restauración a base de Resina </div>	<div style="display: flex; justify-content: space-between;"> Instructions for use Gebrauchsanleitung </div> <div style="display: flex; justify-content: space-between;"> Mode d'emploi Istruzioni d'uso </div> <div style="display: flex; justify-content: space-between;"> Instrucciones de uso </div>
	

ENGLISH

Please read all Information, Precautions, Notes and materials before using BISTITE II DC.

CONTENTS

BISTITE II Paste (A&B)	3.5g each contained in the Super Syringe
PRIMER 1(A&B)	5mL each
PRIMER 2	5mL
AIR BARRIER	6mL
TOKUSO CERAMIC PRIMER (A&B)	1mL each
METALTITE	1mL
Accessaries : Mixing Well, Mixing Pad, Spatula, Disposable brush & handle, AIR BARRIER Dispensing Tips	
Refill:	
BISTITE II Paste (A&B)	3.5g each contained in the Super Syringe
PRIMER 1A	5mL
PRIMER 1B	5mL
PRIMER 2	5mL
AIR BARRIER	6mL (with AIR BARRIER Dispensing Tip x 50)
AIR BARRIER Dispensing Tip	50 tips
TOKUSO CERAMIC PRIMER (A&B)	5mL each
METALTITE	5mL

CLINICAL FEATURES

- BISTITE II DC is a dual-cured adhesive resin cement system with excellent handling and adhesive properties.
- The main components of the BISTITE II DC kit are the resin cement pastes contained in the Super Syringe and a series of primers used to enhance bonding to various substrates.
- Super Syringe allows easy handling by extruding base and catalyst equally at the same time.
- The primer components in the kit consist of PRIMER1 and PRIMER2 for teeth bonding.
- METALTITE, which contains a thiouracil monomer as an adhesive promoter, enhances the adhesion of resins to precious metal alloys with no tin plating.
- TOKUSO CERAMIC PRIMER , which contains a silane coupling agent, enhances bonding between ceramics and resins without hydrofluoric acid etching.

<p>CHEMICAL COMPONENTS</p>	<p>Components</p>
Resin Cement Pastes (A and B)	Filler: Silica-zirconia (77wt%)
PRIMER 1 (A and B)	Matrix: Dimethacrylate, MAC-10(an adhesive promoter), initiator
PRIMER 2	Phosphoric acid monomer, acetone, alcohol, water, initiator
PRIMER 2	HEMA, acetone, initiator
METALTITE	MTU-6 (an adhesive promoter), alcohol
TOKUSO CERAMIC PRIMER (A and B)	Silane coupling agent, phosphoric acid monomer, alcohol

CLINICAL APPLICATIONS

- Cementation and bonding of precious metal, non-precious metal, cured composite resin, porcelain and ceramic prosthesis.
- Repair of fractured porcelains fused to metal crowns and all ceramic
- Cementation of acid etch bridges.
- Cementation of resin composite core post systems and metal post and cores.
- NOTE: It is important to perform occlusal adjustment to avoid fracture by excess occlusal pressure, especially in the case of ceramic occlusal surfaces on posterior teeth.

CONTRAINDICATIONS

*BISTITE II DC contains methacrylic monomers. BISTITE II DC is contraindicated in patients allergic to or hypersensitive to methacrylic and related monomers and should NOT be used in such patients.

CLINICAL PROCEDURES

1 . Tooth Preparation

- Remove contaminants on the tooth surface such as plaque by brushing with flour of pumice or ultrasonic scaling.
- Further, clean the surface of the tooth with alcohol, citric acid, EDTA, or phosphoric acid to improve adhesion and to remove any of the following:
 - Temporary cement residue.
 - Oil from materials to test crown fit.
 - Oil mist from a hand-piece.
 - Saliva, blood, exudate fluids.
- Mixing the PRIMER 1A and 1B
 - Place one drop each of PRIMER 1A and 1B onto a same hole on the mixing well.
 - Mix thoroughly and evenly with the mini sponge or the disposable brush.
- Application of PRIMER1
 - Apply the mixed PRIMER1 immediately onto both enamel and dentin cavity surface with the sponge or the brush.
 - Let it sit for 30 seconds on the surface of the tooth being prepared. On vertical cavity walls the PRIMER1 might settle (lower) or drip (upper). To provide proper treatment , apply the PRIMER1 several times during the 30 second application time.
- Several air dry for 2 to 3 seconds. DO NOT RINSE.
- Application of PRIMER2
 - Place one drop of PRIMER 2 onto another hole on the mixing well. DO NOT mix PRIMER 2 with PRIMER 1 in the hole.
 - Apply PRIMER 2 on the treated tooth surface previously coated with PRIMER 1.
 - Let it sit for 20 seconds on the surface.
- Gently air dry for 3 to 5 seconds.

NOTE: Always apply PRIMER 1A and 1B first, then PRIMER2. Never allow PRIMER 1A and 1B or PRIMER2 to stay over three minutes in the mixing well.

Evaporation may diminish their effectiveness.

I . Restoration Pretreatment

Follow the each procedure shown below, depending on to the materials to be restored:

- Non-precious Metal Prosthesis**
 - Sandblast the surface to be cemented of prosthesis thoroughly with alumina powder (30 to 50 µ m).
 - Ultrasonically clean the restoration.
 - Allow it to dry.
- Ceramics**

*Be sure to read PRECAUTIONS 1.

 - Roughen the surface of the prosthesis with sandblast or grinding with a diamond bur using a slow speed to prepare the area for adhesion.
 - Ultrasonically clean the prosthesis.
 - Allow it to dry.
 - Place a drop of each TOKUSO CERAMIC PRIMER A and B onto a same hole on the mixing well.
 - Mix thoroughly and evenly with the disposable brush or the mini sponge.
 - Apply the mixed TOKUSO CERAMIC PRIMER on the ceramic surface to be bonded thoroughly. DO NOT over apply.
 - Allow the material to stand for 10 seconds prior to bonding.

- NOTE: This procedure is when the ceramic is not cemented on the abutment tooth.
- If the ceramic is on the tooth, clean the prepared surface with non-oily flour of pumice, rinse with water, dry then apply TOKUSO CERAMIC PRIMER.
- Proceed with the bonding procedure according to the manufacturer’s instruction for the use of the adhesive material.
- Resins**
 - Roughen the surface of the prosthesis with sandblast or a diamond bur using a slow speed to prepare the area for adhesion.
 - Follow the same procedure used for ceramics.
 - Precious Metal Prosthesis**
 - Roughen the surface of the prosthesis with sandblast or grinding with a green stone. When using the sandblast procedure, use at least 30 to 50 micron grit.
 - Ultrasonically clean the restoration.
 - Allow it to dry.
 - Apply METALTITE on the surface to be bonded. One or two applications is sufficient. Too thick a layer may result in poor adhesion. Blow drying is not necessary. The solvent of METALTITE evaporates rapidly.

NOTE: DO NOT let saliva, blood, grease, water etc., come in contact with the surface to be bonded. If the surface is contaminated, wipe the surface with an alcohol soaked mini sponge or cotton swab to remove contaminants.

There is no need to reapply METALTITE as it cannot be removed with alcohol.
 - Bond the surface with BISTITE II Paste. Cover the surface as quickly as possible, after the surface has been covered with METALTITE.

II . Cementation and Final Bonding

- How to use the Super Syringe**
 - Lift the cap with the lever.
 - Rotate the lever fully. Always use the lever to lift or rotate the cap.
 - Place the nozzles onto the mixing pad.
 - Push the center of the plunger to extrude the same amount of each paste.
 - Wipe and clean the tips of the syringe nozzles.
 - Return the lever to the proper position where the nozzles are completely closed. Close the cap to the point you can hear a “click”.
- Mixing and Applying Paste A and B**
 - Use Super Syringe to extrude adequate amount of Pastes A and B onto the mixing pad, and mix evenly for 10 seconds.
 - First, apply the mixed cement on the surface of the prosthesis. Place the prosthesis.
 - Wipe the tips of the syringe nozzles and close the cap of the Super Syringe properly after use.
 - For successful cementation, remind quick operation within the working time.
 - For cast post and core, do not apply the cement in root canals. Apply the cement only on the metal core post and sit it on the root canal with gentle pumping action, then sit the post and core in place.
 - For cementing a ling bridge with multiple abutments, apply the cement to the larger crowns first, then the smaller crowns for proper setting time.
 - The cement cures faster when not exposed to a rich oxygen atmosphere.
- Cementing the Restoration for Proper Setting Time**
 - Place the restoration on the tooth with a pressure.
 - Use excess cement, light irradiate for 20 sec. along marginal lines.
 - In case light irradiation is difficult to perform, apply AIR BARRIER with the Dispensing Tip around the margin of the restoration.
- Final Adjustment**
 - Allow the cement to cure for 3 minutes with the patient applying normal occlusal pressure on the restoration.
 - In the case where the AIR BARRIER is applied, remove the AIR BARRIER and rinse with water after 3 min. of application.
 - Check patient’s bite, if necessary.
 - Repolish the restoration with material of choice.

TECHINICAL INFORMATION

Adhesive Strengths (MPa)				
	Bovine enamel (primed)	Bovine dentin (primed)	Co-Cr Alloy (sandblasted)	An-Pd Alloy (sandblasted and primed with METALTITE)
BISTITE II DC	23.0	15.3	22.8	22.9

Working Time	23°C	37°C	
	4 min	-	
Setting Time	8 min 40 sec	3 min	
Physical Properties			
Compressive Strength (MPa)	Hardness (HK)	Tooth brush abrasion (µ m)	Film Thickness (µ m)
	348	50	2

- PRECAUTIONS**
- When using BISTITE II DC for cementing ceramic restorations, be sure to confirm whether or not the restorations are made of Alumina and/or Zirconia. DO NOT APPLY TOKUSO CERAMIC PRIMER for Alumina and/or Zirconia restorations. It will adversely affect optimum cementation of these restorations. For Alumina and/or zirconia restorations, simply roughen restoration surfaces by using a sandblaster or diamond bur before cementing.
 - Use BISTITE II DC only as directed. DO NOT use BISTITE II DC for any other application.
 - BISTITE II DC is intended for sale and use only by licensed dental care professionals.
 - DO NOT use this product if the safety seals are broken or if the contents have been tampered with.
 - DO NOT use BISTITE II DC directly on the cavity if the cavity has pulpal exposure. Protect pulp with a lining or a base material, if the cavity has pulpal exposure.
 - Use examination gloves at all times while using BISTITE II DC to avoid the possibility of allergic reactions from methacrylic monomers.
 - If patient experiences an allergic reaction or sensitivity to BISTITE II DC, such as dermatitis eruption, discontinue its use immediately.
 - AVOID contact of BISTITE II pastes, BISTITE II primers, METALTITE and TOKUSO CERAMIC PRIMER with mucosal membrane, eyes, skin and clothing.
 - AVOID CONTACT WITH EYES: If BISTITE II pastes, BISTITE II primers, METALTITE and TOKUSO CERAMIC PRIMER come in contact with eyes, thoroughly flush eyes with water and immediately contact a physician.
 - If BISTITE II primers comes in contact with mucosal membrane, wipe it away immediately from the effected area after placement of prosthesis. Sensitive areas may whiten from protein coagulation, but such whitening should disappear within 24 to 48 hours.
 - If BISTITE II pastes, BISTITE II primers, METALTITE and TOKUSO CERAMIC PRIMER comes in contact with clothing or skin even when covered with clothing, immediately wipe the area with alcohol soaked cotton swab or gauze.
 - To avoid cross infection, DO NOT reuse disposable brushes and AIR BARRIER Dispensing Tips. Clean the mixing well and spatula with alcohol after use.

- DO NOT leave BISTITE II DC within the reach of patients and children.
- Open the primer bottles carefully. The internal vapor pressure might splash outward.
- Be sure to screw the right cap to the right bottle after each use.
- DO NOT remove the covers from the BISTITE II primer bottles. They are designed to make placement of primers precise and smooth by protecting the direct heat transfer to bottle.
- DO NOT mix PRIMER 1 and PRIMER 2 in the mixing well.
- DO NOT use MIXING well supplied by other manufacturers. The mixing well could be affected by BISTITE II primers.
- DO NOT mix BISTITE II Pastes or BISTITE II primers with other brands of pastes or primers.
- Be sure to set the new tip for AIR BARRIER to its syringe firmly before use and replace with the cap of syringe after use.
- AVOID open flames when using or disposing primers.
- Be sure to use each primer for the proper purpose with the proper procedure.
- DO NOT apply unnecessary force to Super Syringe, especially immediately after taking out of a refrigerator.
- DO NOT use eugenol materials with BISTITE II DC.
- When using a light-curing unit, protective eye shields, glasses or goggles should be worn.
- Be sure to attach Dispensing Tip to AIR BARRIER tightly.

STORAGE

- Store under refrigeration at 0-10°C (32-50°F).
- AVOID direct exposure of BISTITE II pastes to light and heat.
- DO NOT use BISTITE II pastes, PRIMER1, PRIMER2, METALTITE and TOKUSO CERAMIC PRIMER after the indicated date of expiration.

DISPOSAL

- To dispose of BISTITE II pastes, extrude unused pastes from Super Syringe and mix to harden before disposal.
 - The remaining liquid of PRIMER1, PRIMER2, METALTITE and TOKUSO CERAMIC PRIMER should be absorbed in an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.
- *The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure the product is suitable for an appropriate application before use. Specifications are subjected to change without notification.

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