# SILICONE REMOVER

## ENGLISH

READ ALL THE INSTRUCTIONS BEFORE USING SILICONE REMOVER.

### INSTRUCTIONS FOR USE

SILICONE REMOVER is an auxiliary agent used to help to remove silicone-based resilient denture liner from the acrylic denture and to remove silicone rubber impression material from the impression resin tray. SILICONE REMOVER does not inhibit the curing of silicone-based material in relin material and the reimpression.

### PRECAUTIONS

- SILICONE REMOVER should only be used by licensed dental care professionals; dentists, dental hygienists or dental technicians.
- Use SILICONE REMOVER only as instructed herein. Do not use SILICONE REMOVER for any application other than those noted in these instructions.
- Carefully read and follow all manufacturer's instructions prior to using this product.

### Contents

Silicone Remover (7 mL)

### Clinical Application

- Removal of silicone-based resilient denture liner, such as SOFRELINER TOUGH, from dentures.
- Removal of silicone rubber impression materials from the impression resin trays.

### Safety Precautions

- The use of examination gloves (plastic, vinyl, latex) and protective glasses is recommended at all times to avoid the possibility of allergic reactions from acrylic-based monomers. If SILICONE REMOVER comes in contact with the examination gloves, immediately dispose of the examination gloves and wash hands thoroughly with water. In case of contact with skin, wash the skin with soap, and then thoroughly flush with water.
- Ensure that SILICONE REMOVER does not come in contact with patient. In case of contact with the patient's skin, wash the skin with soap, and then thoroughly flush with water.
- AVOID CONTACT WITH EYES. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician for assistance.
- Do not apply SILICONE REMOVER to the patient's mouth or gums. SILICONE REMOVER must be applied only to the denture or impression resin tray and always with the denture or impression resin tray outside of the patient's mouth.
- SILICONE REMOVER is a highly flammable and volatile liquid. Do not use near an open flame. Use it in a well-ventilated area.
- To avoid the risk of cross-infection, do not contact the nozzle of bottle with the denture or relining material.

### Clinical Precautions

- SILICONE REMOVER may spatter during the peeling procedure. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician for assistance. In case of contact with the patient's skin, wash the skin with soap, and then thoroughly flush with water.
- Always use a clean, new plastic cup and brush when using SILICONE REMOVER. Residual dental materials in previously used cups and brushes, such as self curing resin, contain residual catalysts, and may cure SILICONE REMOVER.
- Do not use this product if the safety seal is broken or if the contents have been tampered with.
- Do not immerse the denture or resin tray into SILICONE REMOVER liquid. Otherwise, the surface of the acrylic resin will be unusable.
- Do not use the liquid of self-curing resin instead of SILICONE REMOVER. Additives in the liquid may influence the adhesion to the denture on relining.
- SILICONE REMOVER may not be used for the removal of acrylic resin based soft denture liner and tissue conditioner.

### Storage Instructions


- Store the material at 0-25°C (32-77°F) at all times.
- DO NOT use the material past its expiration date.
- Replace caps tightly when finished.
- To avoid accidental ingestion or spillage, ALWAYS keep out of the reach of children and patients.

### Disposal Instructions

Any remaining SILICONE REMOVER should be absorbed with gauze or cotton, and disposed of in accordance with Federal and/or local regulations. Avoid any inhalation of the volatile substance.

### Directions For Use

- Expose the acrylic surface by cutting the silicon material 2-5 mm in width and 10-20 mm in length with a surgical knife (see Fig. 1).
- Put a few drops of SILICONE REMOVER in the formed concave portion and allow it to sit for 5 - 10 seconds (see Fig. 2).
- Add a few more drops of SILICONE REMOVER, put the instrument or spatula in the boundary between silicone and acrylic resin and then remove the silicone material from the acrylic resin with the instrument or spatula (see Fig. 3). If necessary, add a few more drops of SILICONE REMOVER to the formed concave portion to aid in removal.

### Relining and Reimpression:

Relining or reimpression is performed as follows:

- Remove the surface of adhesion area of the denture or resin tray by using a grinder.
- Allow it to rinse and dry.

The manufacture is not responsible for damage or injury caused by improper use of this product. It is the responsibility of the user to ensure the product is suitable for its intended application.

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**TOKUYAMA DENTAL CORPORATION**

**MATERIAL SAFETY DATA SHEET**

(According to 91/155/EEC and OSHA regulation 29 CFR 1910.1200(g))

**1. Identification of substance**

Product name	SILICONE REMOVER
Description of Product	SILICONE REMOVER is an auxiliary agent used to remove the silicone-based materials from the acrylic resin. This product is classified as non medical device because of the limited application with outside of the patient's mouth even though it is used in a dental clinic or a dental laboratory.
Manufacturer	Tokuyama Dental Corporation 38-9, Taitou 1-chome, Taitou-ku, Tokyo, Japan
EC Representative (Emergency Contact in EU)	Tokuyama Europe GmbH Oststrasse 10, 40211 Düsseldorf, Germany TEL: +49 211 1754480 FAX: +49 211 357379
North America	Tokuyama America, Inc. California TEL: (650) 571 8872 TEL: (877) 378 3548 (Dental Toll-Free)

**2. Composition/Date on components**

Ingredients	%	Cas.No.	Symbols	R-phrases	OSHA PEL
Methyl Methacrylate	>95	80-62-6	F,Xi	R:11-37/38-43	410mg/m <sup>3</sup>

**3. Hazard identification**

Highly Flammable.  
Repeated prolonged contact may cause skin inflammation.

**4. First aid measures**

After skin contact	Wash with soap and plenty of water.
After eye contact	Wash immediately with plenty of water. Get medical attention.
After swallowing	Induce drinking 250mL of water and get medical attention.

**5. Fire fighting measures**

Suitable extinguishing media	Chemical form, carbon dioxide or dry chemical.
Media not to be used	None
Combustion products and resulting gases	Irritant gas may be released during heating or in case of fire.
Special precaution for firefighters	Use protective clothing and a self-contained breathing apparatus.

**6. Accidental release measures**

Personal protection	Not required.
Environmental protection	Do not empty into drains.
Methods for cleaning up	Clean up with an inert absorbent and dispose of it at special waste collection point.

**7. Handling and storage**

Handling	Use in a well ventilated place.
Storage	Store in a cool, dark place(0-25 ).

**8. Exposure controls/personal protection**

In general	Not required.
Respiration protection	Not required.
Hand protection	Wear suitable protective gloves.
Eye protection	Wear suitable protective glasses.
Skin precaution	Wear suitable protective clothing.

**9. Physical and chemical properties**

Appearance	Clear liquid	Explosive properties	1.7-12.5%
Odor	Characteristic Odor	Oxidizing properties	N/A
pH	N/A	Vapor pressure	N/A
Boiling point	100	Relative density	N/A
Melting point	N/A	Solubility(Water)	N/A
Flash point	14.5	Solubility(Solvent)	N/A
Flammability	Highly Flammable	Partition coefficient	N/A
Autoflammability	421	Others	—

**10. Stability and reactivity**

Conditions to avoid	Heat, sparks and open flames.
Materials to avoid	Strong oxidizing agents.
Hazardous reactions	Incomplete thermal decomposition may cause toxic gases.

**11. Toxicological information**

Inhalation	Irritation of respiratory organs.
ingestion	Slight irritation of digestive organs.
Skin contact	Slight irritation, defatting action, repeated prolonged contact may cause inflammation.
Eye contact	Irritation, Redness, Pain.
Carcinogenicity	N/R
Others	—

**12. Ecological information**

Biodegradability	N/A
Bioaccumulative	N/A
Aquatic toxicity	N/A
Others	—

**13. Disposal consideration**

Disposal must be made according to federal, state and local regulations.

**14. Transport information**

Proper shipping name	Methyl methacrylate monomer, stabilized.
UN number	UN1247
Classification	3(IATA,CFR49)
Packing group	II (IATA,CFR49)
Air transport	Highly Flammable (see section 2 and 9). Non medical device (see section 1)
Sea transport	
Inland transport	

**15. Regulatory information**

Marking according to EU regulations	62/548/EEC, 92/32/EEC
Contains	Methyl Methacrylate
Symbols	F: flammable Xi: Irritant
R-phrases	R11: Highly flammable R37/38: Irritating to respiratory system and skin. R43: May cause sensitization by skin contact.
S-phrases	S2: Keep out of the reach of children. S24: Avoid contact with skin. S37: Wear suitable gloves. S46: If swallowed, seek medical advice immediately and show this container or label.
Others	R and S-phrases are not marked according to contents not more than 125mL and the case not related to the general public(Article 23"Labeling"-3).

**16. Other information**

The information which is contained in this document is based on available data. However, as such has been obtained from various sources, including independent laboratories, it is given without warranty or representation that is complete, accurate and can be relied upon. Tokuyama Dental Corp. has not attempted to conceal in any way the deleterious aspects of the product listed herein, but makes no warranty as to such.

\* N/E: Not Established, N/A: Not Available, N/R: Not Regulated by IARC

Manufacturer Tokuyama Dental Corporation 38-9, Taitou 1-chome, Taitou-ku, Tokyo, Japan Tel: +81-3-3835-7201 Email: export@tokuyama-dental.co.jp	EC Rep. Tokuyama Europe GmbH Oststrasse 10, 40211 Düsseldorf, Germany Tel: +49 211 1754480 Fax: +49 211 357379	North America Tokuyama America, Inc. California Tel: (650) 571 8872 Tel: (877) 378 3548 (Dental Toll-Free)
--	---	---