SOFRELINER TOUGH

Silicone-based soft denture liner

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ENGLISH

Carefully read all of these instructions before using SOFRELINER TOUGH.

■PRODUCT DESCRIPTION AND GENERAL INFORMATION

SOFRELINER TOUGH is an addition-cured silicone chairside soft lining material for dentures. SOFRELINER TOUGH has no odor or taste. It is provided in an automix cartridge delivery system, along with an adhesive primer and other accessories. The material cures in 5 minutes intraorally or in 20 minutes at room temperature (23°C/73.4°F).

■CONTENTS:

Kit: PASTE / Base (27g) and Catalyst (27g)×1
Primer (10mL)×1
Drip-Cap for Primer×1
TD II Mixing Tip (XS)×10
Spatula (No. 001)×1
Brush (No.5)×2
Plastic Cup×2
Marking Pen×1
SHAPE ADJUSTMENT POINT (White)×1
FINISHING POINT (Brown)×1
DIAMOND ROUND BUR (for relief)×1
* Note: Use a dispenser that fits the cartridge of SOFRELINER TOUGH.

Refill: Paste / Base (27g) and Catalyst (27g)
Primer (10mL)
TD II Mixing Tip XS
TD II Mixing Tip S
SHAPE ADJUSTMENT POINT (White)
FINISHING POINT (Brown)
DIAMOND ROUND BUR
* Please check with local distributor for refills

■CLINICAL APPLICATIONS FOR SOFRELINER TOUGH

- For pain caused by sharp ridges of alveolar bone and mucous atrophy.
- Poor retention caused by extreme ridge atrophy
- Major undercuts or maxillary/mandibular tori
- Retention of overdentures before the final attachment is placed

■CONTRAINDICATIONS

SOFRELINER TOUGH is contraindicated in patients allergic to or hypersensitive to silicone and platinum compounds and should NOT be used on such patients.

■PRECAUTIONS

- SOFRELINER TOUGH is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- Use SOFRELINER TOUGH only as instructed herein. DO NOT use SOFRELINER TOUGH for any application other than those noted in these instructions.
- DO NOT use SOFRELINER TOUGH if the safety seal is broken or if it appears the contents have been tampered with.
- Primer contains ethylacetate, which is a highly flammable and volatile solvent. Do not expose the primer to open flame. Use in a ventilated room during use and close the cap tightly after use.
- If patient experiences an allergic reaction or sensitivity to SOFRELINER TOUGH, such as dermatitis eruption, remove the prosthesis and discontinue its use immediately.
- Wear non-latex examination gloves at all times while using SOFRELINER TOUGH to avoid the possibility of allergic reactions. Latex gloves may also inhibit the cure of SOFRELINER TOUGH. If SOFRELINER TOUGH primer comes in contact with the examination gloves, immediately dispose of the examination gloves and wash hands thoroughly with water.
- Avoid contact with eyes. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician.
- Avoid contact of uncured chemical with skin. In case of contact with skin, clean affected area with rubbing alcohol, then thoroughly flush with water.
- Avoid contact of unpolymerized materials with clothing. In case of contact with clothing, do not wipe the material. Allow it to set and then peel it off.
- DO NOT allow excess materials to flow down the patient's trachea or esophagus. Aspirated or ingested SOFRELINER TOUGH paste can cause asphyxiation or other injury. Immediately remove any excess lining material from the palatal area to avoid material flowing down the patient's throat. SOFRELINER TOUGH is not radiopaque, and x-rays will not detect aspirated or ingested material.
- DO NOT mix SOFRELINER TOUGH with any materials not noted in these instructions.
- After using SOFRELINER TOUGH, the mixing tip must be replaced with the cartridge cap that comes with the cartridge at the time of purchase. Clean the tip of the cartridge with a clean cloth or tissue each time before using a new mixing tip or replacing the cartridge cap.
- The material cures in 5 minutes intraorally or in 20 minutes at room temperature (25°C / 77°F). SOFRELINER TOUGH will be completely cured in 24 hours after the lining. Note that the softness of the material may change slightly during this period. The finishing procedure can be performed immediately after the extraction of the newly lined denture from the mouth or a gypsum model.
- Always clean the spatula, brushes and cups with alcohol after each use to avoid cross contamination.
- Instruct patients in accordance with the PATIENT PRECAUTIONS contained in the last section of these instructions.

- To avoid inhaling dust while polishing, always wear a mask recommended by federal and/or local regulations.
- Read and follow the manufacturers' instructions for all products mentioned in this instruction manual that are not included in this kit.
- As much as possible, avoid contact with oral mucous membrane other than the region to be relined. Instruct patients to rinse their mouth with adequate amount of water after the procedure.
- When used on patients with reduced reflecting functions because of age, instruct them to breathe nasally and confirm closed respiratory tract.
- Do not bring the Primer into contact with plastic material because it may dissolve some types of plastic.

■NOTES FOR USE

- If there is a great amount of ridge resorption, line the old worn denture first with a new acrylic material, such as self-curing acrylic resin or hard denture liner etc., to enhance the bond between the denture and SOFRELINER TOUGH.
- Use SOFRELINER TOUGH Primer to bond acrylic surface. Use only the Primer supplied with SOFRELINER TOUGH, and use only a thin coat of Primer. DO NOT apply SOFRELINER TOUGH Primer to the patient's mouth or gums. SOFRELINER TOUGH Primer should be applied only to the denture and always with the denture outside of the patient's mouth.
- If the SOFRELINER TOUGH Primer may become thickened because of volatilization, do not use it as the adhesive strength may be reduced.
- Always use a mixing tip to dispense. Do not mix pastes by hand because air bubbles may be incorporated into the paste, which may cause deterioration, stain and unpleasant odor.
- To achieve a thicker surface, allow the SOFRELINER TOUGH Paste dispensed onto the denture to stiffen slightly before inserting the denture into the patient's mouth.
- Avoid contact with the materials below, as they may inhibit the curing of SOFRELINER TOUGH:
Silicone-based impression materials, eugenol containing materials, polysulfide rubber impression materials, natural rubber gloves, other manufacturer's soft lining materials and uncured acrylic resins.
- Soaking the denture in heavy bleach cleaning solutions may alter the color of SOFRELINER TOUGH, but it does not affect the use of denture.
- Instruct patient to clean the denture as described in PATIENT PRECAUTIONS.
- Since silicone-based fitting test material sticks to SOFRELINER TOUGH, apply vaseline on the surface prior to use.
- If there is a fracture in finished reline, remove all the relining completely and then perform a new procedure.
- Refer to following materials and methods for shape adjustment and polishing.

Stages	Materials	Shapes
Trimming	Sharp scalpel ¹⁾	Slightly rounded blade ²⁾
Shape Adjustment	SHAPE ADJUSTMENT POINT (White) ²⁾	Round end taper
Finishing	FINISHING POINT (Brown) ³⁾	Round end taper
Inner Surface Adjustment	DIAMOND ROUND BUR ⁴⁾	Round

- Notes : 1) Strong trimming by scalpel is easy to make sharp notches that cause fracture.
The sharp notches should be made smooth by polishing.
2) Correct the form by sliding the point lightly on the surface under 15,000 rpm.
3) Finish by sliding the point lightly on the surface under 4,000 rpm.
4) Use DIAMOND ROUND BUR under 15,000 rpm.
5) Surgical blade No. 10.

- Excessive removal of denture resin could cause denture fracture and weakening.

■CLINICAL APPLICATION OF SOFRELINER TOUGH INTRAORAL LINING METHOD


To achieve the maximum effect using SOFRELINER TOUGH, the optimum thickness of lining material should be 1-2 millimeter. Areas which will require lining thicker than 2 millimeters should first be lined with TOKUYAMA REBASE II, chairside hard denture lining material, to obtain the proper thickness for SOFRELINER TOUGH. Material thickness of SOFRELINER TOUGH will not affect its curing.

1. EXAMINATION OF MOUTH

Check intraoral condition by inspection and palpation.


2. LINING

- Remove denture from patient's mouth and clean the area to be lined with soap and water. Rinse thoroughly, and then air-dry the denture surface thoroughly.
- Mark a red borderline on the denture with the Marking Pen (See Figure 1) in the case that the paste is to be placed on the labial and/or buccal side passing over the denture border. Avoid setting the marginal line on the denture border.


- Grind the tissue side of the denture with a bur to the proper depth. Bevel the borders of the denture to allow 1 millimeter of material to roll over the borders (See Figure 2).

Note: When the border is made within the flange of denture base, remove resin in the relining region by using round bur (for acrylic resin) to ensure proper resin depth.

- Blow the ground denture surface with strong air stream. Pour the Primer into the Plastic Cup. Using the Brush, apply a thin, even coat of Primer on the denture surface where SOFRELINER TOUGH is to be lined. Do not miss any spot where the Primer should be applied since this may cause poor adhesion. Blow air to dry Primer, or leave it in well-ventilated area for 20 seconds.
Note: Attach the Drip-Cap onto the rim of the Primer bottle to facilitate pouring. Primer application on the red borderline will blur the color slightly. To avoid blurring the red borderline, try to minimize the amount of primer applied.
- Load the SOFRELINER TOUGH paste cartridge to a dispenser (sold separately). Attach Mixing Tip to the cartridge by aligning the notches. Then, push it onto the cartridge and then rotate the Mixing Tip 90 degree clockwise.
- Dispense SOFRELINER TOUGH paste evenly on the tissue side of the denture. Using a spatula, spread the paste evenly across all surfaces where the Primer was applied.
Note: If any part of the primed surface is not covered with paste in this stage, the area may be contaminated by saliva when it is inserted into the mouth and thus cause low adhesion to the paste in that area.
- Suction any excess saliva from patient's mouth within 1 - 1 1/2 minutes after dispensing SOFRELINER TOUGH onto the denture. Place the denture in patient's mouth, and have the patient occlude very gently into centric occlusion with the proper vertical dimension. Check for any excess SOFRELINER TOUGH in the patient's mouth, and remove immediately if present. DO NOT allow excess material to flow down the patient's trachea or esophagus. Aspirated or ingested SOFRELINER TOUGH paste can cause asphyxiation or other injury, even though SOFRELINER TOUGH paste is non-toxic. If the material flows towards the patient's trachea, immediately wipe excess material from patient's mouth with a cotton ball or gauze. Advise patients to alert the dentist any time he/she feels SOFRELINER TOUGH material or any substance flowing toward the throat.
Note: To line upper dentures with SOFRELINER TOUGH, set the denture from the posterior palatal region to avoid flow of SOFRELINER TOUGH material down patient's throat.
Take care not to place too much SOFRELINER TOUGH paste on the denture and place it earlier in the case of relining the palatal region because hardened paste tends to become thick in this case.

- Pull the cheeks very gently to border mold the vestibule area while the patient is gently occluding.
 - Allow the denture to completely set in patient's mouth while monitoring the patient for any excess SOFRELINER TOUGH. If patient salivates, suction the saliva and instruct to resume occluding gently.
 - Since silicone-based fitting test material sticks to SOFRELINER TOUGH, apply vaseline on the surface prior to use.
 - If there is a fracture in finished reline, remove all the relining completely and then perform a new procedure.
 - Refer to following materials and methods for shape adjustment and polishing.
- Note: DO NOT over trim the area where the bevel was made on the denture. DO NOT leave any loose SOFRELINER TOUGH material on the denture.

3. FINISHING

Trim excess with sharp scalpel as much as possible and use SHAPE ADJUSTMENT POINT (White) to make the surface smooth or to adjust the shape.
When necessary, use the FINISHING POINT (Brown) for finishing (See Figure 4).

Note: FINISHING POINT (Brown) should be used by rotating from relining area to denture area in order not to break the margin region (See Figure 4).
Use DIAMOND ROUND BUR or any polishing material suitable to make a relief treatment. Erase the red marginal line with FINISHING POINT (Brown) or others.

4. OTHER IMPORTANT CONSIDERATIONS FOR THE SUCCESSFUL RELINE

In order to enhance bonding to old or worn dentures, line the denture first using TOKUYAMA REBASE II. Allow it to cure completely prior to lining with SOFRELINER TOUGH. Make sure that the tissue side of the denture is ground deeply enough, so that the relining of SOFRELINER TOUGH will not be too high on the patients' occlusion after the lining with TOKUYAMA REBASE II.

■CLINICAL APPLICATION OF SOFRELINER TOUGH LABORATORY LINING METHOD

1. EXAMINATION OF MOUTH

Check intraoral condition by inspection and palpation. Examine the area where needs relief treatment.

2. EXAMINATION OF THE DENTURE

If the denture is worn or old, line the denture first using TOKUYAMA REBASE II, chairside hard denture lining material. The hard lining material should be applied in the appropriate thickness so that when the new material is ground down, a newly lined layer of hard lining material is left on the denture bearing area, allowing SOFRELINER TOUGH to adhere to it. Make sure that the tissue side of denture is ground deeply enough, so that the relining of SOFRELINER TOUGH will not be too high on the patients' occlusion after the lining with TOKUYAMA REBASE II.

3. FABRICATION OF A GYPSUM MODEL

Take an impression with the patient biting on the denture, then pour a stone model on the impression.

- Make a stone index of the teeth and mount the denture on the

- index of a lining jig.
- Make a relief procedure in the stone model with wax in the reign where may cause pain in the mouth.
- Mark a red borderline denoting the area to be relined on the denture with the Marking Pen (See Figure 1).
- To ensure proper thickness of the lining material, grind the lining area with an acrylic bur (See Figure 2).
Note: Use of the round bur makes it easy to grind by measuring the depth with the diameter of the bur as the guide.
- Replace the denture on the index and position the denture on the lining jig.
- Apply a separating medium on the ridge area of the gypsum model so the lining material will not adhere to it.
- Apply Primer to the entire area of the denture to be lined.
- Load SOFRELINER TOUGH paste cartridge to a dispenser and apply the paste using Mixing Tip onto the surface where the primer has been applied. Place the paste both on the mucosal side of denture and stone model.
- Place the denture onto its gypsum stone model, then press it with the lining jig.
- Allow denture to sit for 20 minutes or more at room temperature (23°C / 73.4°F) on the lining jig.
Note: Soak the denture in warm water at 40°C to 50°C (104°F to 122°F) for 5 minutes together with the lining jig in order to speed up the curing process.
- After SOFRELINER TOUGH has been cured, remove the denture from the gypsum model, and carefully trim any excess SOFRELINER TOUGH material from the border areas with a sharp scalpel (See Figure 3). DO NOT over trim the area where the bevel was made on the denture. DO NOT leave any loose SOFRELINER TOUGH material on the denture.

4. FINISHING

Trim excess with sharp scalpel as much as possible and use SHAPE ADJUSTMENT POINT (White) to make the surface smooth or to adjust the shape (See Figure 4).
When necessary, use the FINISHING POINT (Brown) for finishing (See Figure 4).

■PATIENT PRECAUTIONS

Patients should be informed of the following instructions to enhance the longevity, durability and safety of a denture lined with SOFRELINER TOUGH material:

- SOFRELINER TOUGH lining material should not be scrubbed with a denture brush or toothbrush. The best way to clean SOFRELINER TOUGH is with sponge or a soft piece of gauze. While running cold water on the SOFRELINER TOUGH surface, wipe the surface with sponge or gauze saturated with soap. Use a gentle wiping motion to remove all debris and food deposits from the SOFRELINER TOUGH lining material. The tooth surface and hard acrylic denture base may be cleaned with a regular denture brush by scrubbing with soap, denture toothpaste or regular denture cleaners at the dentist's discretion.
- Soak the denture at least once a day, preferably at bedtime, in an oxygenating denture cleaner (i.e. Polident, Efferdent, etc.). The denture should not be cleaned with an oxygenating cleaner within the patient's mouth. Color fading may occur with long term use, depending upon the patient's denture maintenance. Use a denture cleaner according to its manufacturer instructions.
- Soaking the denture in heavy bleach cleaning solutions may alter the color of the SOFRELINER TOUGH, but it will not affect its performance.
- The denture should be stored in water when not being worn.
- The lined denture should be inspected and ultrasonically cleaned by a dentist every three months, or whenever plaque deposit is noted.
- The SOFRELINER TOUGH must be inspected frequently for any damage. A denture with disturbed, broken, cracked, loose or damaged liner should not be worn. Use of a damaged or broken denture could result in choking, aspiration, asphyxiation or other injury.
- In the event of any sign of sensitivity to silicone, use of the denture should be discontinued, and the patient should contact his/her dentist immediately.

■STORAGE

- Store SOFRELINER TOUGH at a temperature between 0-25°C (32- 77°F).
- Avoid direct exposure of SOFRELINER TOUGH to heat, open flame, humidity and sunlight.
- Do not use SOFRELINER TOUGH after the indicated date of expiration.
- DO NOT leave SOFRELINER TOUGH unnecessarily within the reach of patients and children to avoid being swallowed by mistake.

■DISPOSAL

- To dispose the paste, dispense unused paste out from the cartridge and harden it before disposal.
- The remaining liquid of the primer should be absorbed in an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

*The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure the product is suitable for an appropriate application before use.
Specifications are subject to change without notification.

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