ONE-UP BOND F Plus

CE
0123
REF AD01

PLEASE READ ALL INFORMATION, PRECAUTIONS AND OTHER MATERIALS ACCOMPANYING THIS PRODUCT

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ENGLISH

■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

1. ONE-UP BOND F Plus is a self-etching dental adhesive that permits a simplified clinical procedure eliminating the need to pretreat a tooth when performing direct composite restorations. It may be applied on exposed cervical dentin and the root surface.
 2. ONE-UP BOND F Plus does not require a separate etching or priming step.
 3. ONE-UP BOND F Plus does not require water rinsing or air blow for drying.
 4. ONE-UP BOND F Plus has built-in color-coding. The color of ONE-UP BOND F Plus will turn from yellow to pink upon mixing, and the pink color will disappear upon light irradiation.
 5. ONE-UP BOND F Plus releases fluoride.

■ CLINICAL APPLICATION AND USAGE

1. Rinse the cavity preparation. In case where no cavity preparation has been made, clean the tooth surface with a rubber cup and a fluoride-free cleaning paste; rinse thoroughly with water.
 2. Remove water by air blowing to keep the cavity in a dry or a moist condition whose surface should not show the visible puddles or any shiny appearances.
 3. Using Dispenser provided with ONE-UP BOND F Plus, place one drop of BONDING AGENT A and BONDING AGENT B into the mixing well.
 NOTE
 • Before dispensing BONDING AGENT B, shake the bottle lightly, because the filler in BONDING AGENT B may sediment on the bottom of the bottle during periods of storage.
 • Hold each bottle of BONDING AGENT vertically to ensure that a full drop of BONDING AGENT is dispensed from the bottle. If each bottle is not held vertically, the amount of one drop will be larger than a normal drop and it may also foul the nozzle.
 • If any BONDING AGENT A and/or BONDING AGENT B remains on the tip of the nozzle, wipe the excess.
 4. Mix thoroughly and evenly until the mixed BONDING AGENT turns completely and evenly pink.
 5. Close Dispenser lid tightly after use.
 6. Apply the mixed BONDING AGENT to the cavity, both the enamel and dentin, with the Disposable Applicator Tip.
 NOTE
 • Use mixed BONDING AGENT within 1.5 minutes after mixing. The mixed BONDING AGENT will thicken, even without light irradiation, due to the self-curing characteristic of the mixed BONDING AGENT after mixing.
 • Avoid rinsing or hard air blowing after applying the mixed BONDING AGENT. If saliva, blood, or other fluid contaminates the area where the mixed BONDING AGENT is applied, thoroughly rinse the cavity, dry and re-apply newly mixed BONDING AGENT.
 7. After applying mixed BONDING AGENT, keep rubbing the mixed BONDING AGENT on the surface for at least 10 seconds. ONE-UP BOND F Plus does not require air blowing because it contains no volatile solvent, e.g. ethanol or acetone.
 NOTE
 • If it is difficult to rub the mixed BONDING AGENT on the bonded surface, allow it for at least 20 seconds after applying mixed BONDING AGENT for BONDING AGENT to adhere. If mixed BONDING AGENT is applied two or three times on

the same surface, it will help to ensure proper bonding.

- If the mixed BONDING AGENT becomes overly concentrated in one area of the cavity preparation, blow air gently to spread BONDING AGENT more evenly.
- Do not rinse the applied layer by water.

- If any part of the applied layer of BONDING AGENT flows and becomes thin, reapply BONDING AGENT. The adhesive strength becomes lower when the applied layer becomes too thin. Maintain the thickness of the applied layer on condition that the glossiness of the surface can be observed (more than 10 µm in thickness).

8. Light irradiate the mixed BONDING AGENT for at least 10 seconds keeping the curing light tip within about 2 mm from the cavity surface.

NOTE

- The pink color should turn to a pale brown after light irradiation.

- Be sure to confirm that the pink color has disappeared upon light irradiation. In case the color of BONDING AGENT remains pink, irradiate with light for an additional period of time until the pink color disappears.

9. Restore with composite resin (do not use chemical-cured composite resins).

NOTE

- DO NOT shape composite resin paste with any filling instrument that is coated with mixed BONDING AGENT. The paste surface may otherwise be stained by the color of mixed BONDING AGENT.

■ PRECAUTIONS

1. ONE-UP Bond F Plus is intended for sale to, and use only by, licensed dental care professionals.
2. Use ONE-UP BOND F Plus only as directed. DO NOT use ONE-UP BOND F Plus for any other applications.
3. The use of examination gloves (plastic, vinyl, or latex) is recommended at all times to avoid the possibility of allergic reactions from methacrylic monomers.

If ONE-UP BOND F Plus comes in contact with the examination gloves, remove and dispose of the gloves and wash hands thoroughly with water as soon as possible.

4. Always determine whether your patient is allergic or hypersensitive to methacrylic monomers, organic solvents, and/or acids before using ONE-UP BOND F Plus. DO NOT use ONE-UP BOND F Plus if the patient is allergic to methacrylic monomers, organic solvents, and/or acids.

5. If ONE-UP BOND F Plus causes an allergic reaction or oversensitivity, discontinue its use immediately.

6. Avoid contact of ONE-UP BOND F Plus with eyes, mucosal membrane, skin and clothing.

- If ONE-UP BOND F Plus comes in contact with eyes, thoroughly flush eyes with water and immediately contact a physician.

- If ONE-UP BOND F Plus comes in contact with mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration.

- Affected areas may whiten from protein coagulation, but such whitening should disappear within 24 hours. If such whitening does not disappear within 24 hours, immediately contact a physician, and the patient should be so advised.

- If ONE-UP BOND F Plus comes in contact with clothing or sensitive skin, even when covered with clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.

- 7. Use only Disposable Applicator Tip provided with ONE-UP BOND F Plus, as these tips are designed to yield the precise amounts of bonding agents required.

8. DO NOT expose the product or its vapor to open flame.

9. Open BONDING AGENT bottles carefully. The internal vapor pressure might cause splashing of BONDING AGENT.

10. DO NOT leave ONE-UP BOND F Plus within reach of patient or children.

11. DO NOT use this product if the safety seals are broken or if the contents appear to have been tampered with prior to use by the dental professional.

12. To avoid cross infection, DO NOT reuse Disposable Applicator Tip.

13. Clean the mixing well thoroughly with alcohol after each use.

14. DO NOT use ONE-UP BOND F Plus with chemical-cured composite resin.

15. DO NOT mix BONDING AGENT with other brands of primer or adhesives.

16. If the cavity is immediately proximate to the pulp, pulp should be protected with a Glass Ionomer lining or calcium hydroxide. DO NOT USE EUGENOL MATERIAL to protect the pulp. Eugenol material will inhibit curing of ONE-UP BOND F Plus.

17. When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.

■ STORAGE

1. Store under refrigeration at a temperature between 0-10°C(32-50°F).
2. DO NOT use ONE-UP BOND F Plus after the indicated date of expiration.

■ DISPOSING OF MATERIALS

Unused BONDING AGENT and constituents A and B should be absorbed into inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

*The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure the product is suitable for an appropriate application before use.

Specifications are subjected to change without notification.

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