

ONE-UP® BOND F Plus

Self-Etching Bonding Agent


One-Up® Bond F Plus is a light-cured, self-etching bonding agent that permits a simplified procedure by eliminating the pretreatment of the tooth when performing direct composite restorations. It requires 3 simple steps; Mix–Apply–Light cure, and each step can be visually confirmed by the built-in color coding, eliminating guesswork. One-Up® Bond F Plus may be applied on exposed cervical dentin and the root surface.

Benefits

- Visual confirmation of each step by built-in color coding
- Single application
- Fluoride releasing
- No etching, rinsing or drying
- New dispensing system

Indications

- Direct light cured composite restorations.
- Application on exposed cervical dentin and the root surface.



- One-Up® Bond F Plus
- Hybrid Layer
- Dentin

Resin-dentin Interface

ONE-UP® BOND F Plus

Self-Etching Bonding Agent

Product Packages Available

One-Up® Bond F Plus Kit


- Bonding Agent A/5mL
- Bonding Agent B/5mL
- Mixing Well
- Applicator Handle
- Disposable Applicator Tip/50

Refill

- Bonding Agent A/5mL
- Bonding Agent B/5mL
- Disposable Applicator Tip & Handle Set/500 tips and 4 handles

ONE-UP® BOND F PLUS PRODUCTS

ITEM#	ITEM DESCRIPTION
14706	One-Up® Bond F Plus Kit
14736	Bonding Agent A
14737	Bonding Agent B
14672	Disposable Applicator Tip & Handle Set



For more information, please call 1.877.ESTELITE or visit www.tokuyama-us.com

