

# PALFIQUE® ESTELITE® LV

Esthetic Flowable Composite Resin

Palfique® Estelite® LV is an esthetic flowable resin composite, which offers 3 different viscosities; High Flow, Medium Flow and Low Flow, to suit different clinical cases and the preferences of the clinicians. Every inorganic filler contained in Palfique® Estelite® LV is a submicron spherical filler (0.4um and 0.08um), providing the same benefits as Estelite® Sigma, such as excellent esthetics, great self-shade matching, superior polishability and gloss retention etc. Palfique® Estelite® LV Opaquer and Palfique® Estelite® LV Clear are also available.




## Benefits


- Offering 3 different flowabilities
- Outstanding esthetics
- Excellent polishability
- Exceptional surface gloss
- Super low shrinkage
- Wider shade matching range / "Chameleon effect"
- Easy "non-sticky" handling
- High wear resistance
- Optimal radiopacity
- Minimum color change after curing

## Indications

- Restoration of class-V cavities
- Restoration of deciduous teeth
- Restoration of incipient caries
- Restoration of small class-I, class-II and class-III cavities
- Restoration of small cavities or tunnel preparation
- Basing before composite resin restoration
- Lining of class-V cavities prior to restoration with other composite resins
- Contour correction after composite resin filling
- Repair of fractured porcelain
- Repair of small fractured composite resin
- Blocking of undercuts




Palfique® Estelite® LV



Filtek™ Flow\*



Tetric® Flow\*



Revolution® Formula 2™\*

\*Not a registered trademark of Tokuyama Dental Corp.

# PALFIQUE® ESTELITE® LV

Esthetic Flowable Composite Resin

## Product Packages Available

### Palfique® Estelite® LV Syringe

- 3g/syringe + 15 syringe tips


### Palfique® Estelite® LV Syringe Tip

- 50 tips/box



## Flowability and Shades

- **High Flow:**  
A1, A2, & A3
- **Medium Flow:**  
A1, A2, A3 & A3.5
- **Low Flow:**  
A1, A2 and A3



## PALFIQUE® ESTELITE® LV PRODUCTS

ITEM #	ITEM DESCRIPTION
12341	Palfique® Estelite® LV High Flow Syringe (A1)
12342	Palfique® Estelite® LV High Flow Syringe (A2)
12343	Palfique® Estelite® LV High Flow Syringe (A3)
12351	Palfique® Estelite® LV Medium Flow Syringe (A1)
12352	Palfique® Estelite® LV Medium Flow Syringe (A2)
12353	Palfique® Estelite® LV Medium Flow Syringe (A3)
12354	Palfique® Estelite® LV Medium Flow Syringe (A3.5)
12346	Palfique® Estelite® LV Low Flow Syringe (A1)
12347	Palfique® Estelite® LV Low Flow Syringe (A2)
12348	Palfique® Estelite® LV Low Flow Syringe (A3)
12335	Palfique® Estelite® LV Clear Syringe
12336	Palfique® Estelite® LV Opaquer Syringe
12371	Palfique® Estelite® LV Syringe Tip

For more information, please call 1.877.ESTELITE or visit [www.tokuyama-us.com](http://www.tokuyama-us.com)

